


CLINDIAG
CLINDIAG SYSTEMS B.V.B.A.



SA-10/20

Chemistry Analyzer

NEW



ISO9001:2008 & ISO13485:2003



Headquarters

Add: Steenberg 66, Pollare/Ninove(9401), Belgium
Tel: 0032 54 250 936
Fax: 0032 54 243 058
Website: <http://www.clindiag.be>
E-mail: marketing@clindiag.be

U.S. Office

Add: 1351 S. Leavitt Ave, # 104 Orange City, Florida 32763 USA
Tel: 001 386 456 1235
Fax: 001 386 456 1237
E-mail: info@clindiagusa.com

India Office

Add: B-1177-78, G.D. Colony Mayur Vihar Phase-III, Delhi-110096
PH: 0091 11 22618779
Fax: 0091 11 22618780
E-mail: contact@clindiagindia.com

Features

Technical Specifications

Model	SA-10	SA-20
Reading Cuvette	Direct reading cuvette	Both through cell and direct reading cuvette
Incubating Positions	20 incubating positions	
	Light source: 6V, 10 W longlife halogen lamp	
Photometric System	Filters: 340/405/492/510/546/578/620nm, 2 more open filter positions, ±2nm pass-hand	
	Wavelength accuracy: ±2nm	
Measuring System	Measuring range: 0-2.500 O.D.	
	Photometric linearity: ±2% from 0 to 2.000 O.D.	
	Photometric accuracy: ±2% from 0 to 2.000 O.D.	
	Drift < 0.005 O.D.	
	Carry over ≤ 1%	
Incubator Temperature	25°C, 30°C, 37°C RT	
Control	Precision: ±0.1°C	
Display	5 inch colorful touch screen	
Printer	Built-in thermal printer	
Connection Port	RS-232 & USB	
Power Supply	AC 110V 60Hz / 220V 50Hz	
Dimensions	34cmX38cmX18cm	
Weight	8.5kg	

CLINDIAG...

The Reference for Quality and Service