


CBS-300 Electrolyte Analyzer (K^+ Na^+ Cl^-)



Features

- Advanced Multi-Biosensor Technology (MBT) permits maintenance-free electrode.
- Close reagent pack is easy to replace.
- Multi-level calibrators and controls secure results accuracy.
- User friendly operation and quick trouble shooting.
- High throughput gives each test result within half a minute.
- Up to 1000 sample results storage and easy to review.

Specifications

Method	Multi-Biosensor Technology			
Sample Size	60-120 μ l			
Speed	\geq 80 tests/hour			
	Performance	Detection Range(mmol/L)	Precision(CV)	Resolution(mmol/L)
Measure scope & Precision	K ⁺	0.50 ~ 20.00	\leq 1.0%	0.01
	Na ⁺	15.0 ~ 200.0	\leq 1.0%	0.1
	Cl ⁻	15.0 ~ 200.0	\leq 1.0%	0.1
Printer	Inner thermal printer			
Expand Function	RS232 & COM Connector, Computer and Barcode reader available			
Ambient Conditions	0-40 $^{\circ}$ C, Relative humidity \leq 85%			
Power	AC 80V ~ 260V, 50-60Hz, \leq 80W			

