

Specification

Dimension (mm)	269, 130, 137
Host Weight	1.5 Kg
Data Storage	8000 tests result
Adapter Power	AC100~240V, 50/60Hz
Data output	LIS/HIS, PC, Printer
Rated Power	36W
Displayer	7 inch high resolution touch screen
Battery working time	> 8 hours
QC code reading	RFID
Repeatability	< 3%

Lamuno X

Multi Scenes Usage

Fluorescence Immunoassay Analyzer



AEHEALTH LIMITED

Add: Unit G25 Waterfront Studios, 1 Dock Road, London, United Kingdom, E161AH

Email: sales@aehealth.uk




Designed by "reddot" award winner team  
"Red-Star" award winner team

## Highlight Feature



## Testing Process



## Application scenarios



Drug Rehabilitation Center



Cardiology



ED (Emergency Department)



Clinical Lab



ICU



Surgery



Outdoor Drug Testing Station



Ambulance

## Tests List

### Cardiac Marker

- NT-proBNP
- CK-MB/cTnl/Myo
- cTnl
- Myo
- D-Dimer
- cTnT
- CK-MB
- sST2

### Tumor Marker

- Ferritin
- AFP
- CEA
- CA15-3
- PSA
- f-PSA
- CA125
- CA19-9

### Hormone

- $\beta$ -HCG
- FSH
- E2
- Tes
- LH
- PRL
- AMH
- Prog

### Thyroid Function

- T3
- FT3
- Cor
- T4
- FT4
- TSH

### Drug Abuse

- CAF
- MCAT
- TRA
- KET
- MDMA
- Multi-Drug 3 in 1 test
- THC
- MET
- MOP
- K2
- COC

### COVID-19

- COVID19 Neutralizing Antibody
- COVID19 Antigen
- COVID19 IgM/IgG

### Inflammation Detection

- Hs-CRP + CRP
- SAA
- PCT
- IL-6

### Others

- IgE
- HbA1c
- S100- $\beta$
- 25-OH-VD

### Renal Function

- NGAL
- MAU

### Infectious

- HCV
- HBsAg

### Digestive Function

- G17
- PGI / PGII